

Purpose

This introductory exercise will familiarize you with a few of the measurements we make in the chemistry laboratory and the level of uncertainty associated with each one.

We use quite a few measuring devices in the chemistry lab, but we are typically measuring two simple properties - **mass** and **volume**. Mass measurements are typically made on **balances**, while volume measurements are made on a wide range of devices like **beakers**, **graduated cylinders**, **transfer pipets**, **burets**, and others.

We will look at three of these today - a mass measurement device and two volume measurement devices.

Measuring mass with the analytical balance

For measuring mass, our laboratory is equipped with **analytical balances**. The analytical balance is a high-quality instrument capable of measuring masses with both high accuracy (correctness) and high precision (reproducibility). The trade-off is that you cannot typically weigh large masses on an analytical balance.



Illustration 1 - A typical analytical balance

The balance has a pan, where samples are placed, and is surrounded by a glass cage to prevent air currents from altering the balance reading. Analytical balances usually have simple controls consisting of an on/off button and a "0/T" or "tare" button, which will reset the balance to zero. The balance should be zeroed **before each use**, since all electronic measuring devices have a tendency to "drift"

over time.

If you are weighing an **object** like a beaker (i.e. not a chemical sample), here's how to use the balance.

1. Make sure the balance is on. If the balance is off, turn it on using the power button and **wait 15 minutes** for the balance to stabilize before using it.
2. With all glass doors **shut**, press the "Tare": or "0/T" button to **zero the balance**. All the numbers on the balance display should change to zero. Some of our balances will beep when zeroed, while others may display "Stable". One of our older balances will simply reset the display to all zeros with no other indication that it has done anything.
3. Open one of the glass doors and **set your object in the center of the pan**. **Shut all glass doors**.
4. When the balance reading stabilizes, **write down all digits** displayed by the balance. Also **write down the units**, typically grams.
5. Remove your object from the pan and make sure all the glass doors are shut before leaving.

If you are weighing multiple objects, then **re-zero the balance** for each new object. If you are weighing the same object repeatedly, then **remove the object and re-zero the balance** for each repeat weighing.

If you are weighing a **powder or liquid**, you must protect the balance pan from the substance you are weighing. You can weigh a powder using a sheet of **weighing paper** or a **plastic weighing boat**. For liquids, you may use a small **flask** (50 mL to 100 mL size) or a **plastic weighing boat**.

Here's how to measure out a specific mass of powder.

1. Make sure the balance is on. If the balance is off, turn it on using the power button and **wait 15 minutes** for the balance to stabilize before using it.
2. Place your weighing paper or weighing boat in the center of the balance pan, then shut all glass doors. Press the "Tare" or "0/T" button to **zero the balance**. Zeroing the balance with the paper or boat in place will make the balance **subtract out the paper or boat weight for you** when you add your powder to the paper or boat.
3. Open one of the doors and **carefully add powder to the paper or boat**. Do **not** spill powder on the balance pan. If you're using a weighing boat, it may be easier to simply remove the boat from the balance, add powder, and put the boat back in the center of the pan after adding powder. Add powder until you get the mass you are trying to measure out. If you overshoot the mass, you may gently remove some powder using a spatula. **Do not press down on the balance pan with the spatula**, since a large force on the pan can damage the balance!
4. Close the glass door, wait for the balance reading to stabilize, then **write down all digits** displayed by the balance. Also **write down the units**, typically grams.
5. Remove your paper or boat with your powder and shut the glass doors before leaving. If you have spilled any powder on the balance, gently brush it off the balance using one of the provided brushes. Some powders can corrode the balance pan and destroy the balance if not brushed off!

Do not assume that all papers or boats have the same mass. To weigh out a powder, always zero the balance with the piece of paper or boat that you are going to put your sample on. If you are using paper, use a different piece of paper for each sample. As with weighing objects, re-zero the balance between measurements.

Here's how to measure the mass of a liquid sample.

1. Make sure the balance is on. If the balance is off, turn it on using the power button and **wait 15 minutes** for the balance to stabilize before using it.
2. Place an empty small flask or weighing boat in the center of the balance pan, then shut all glass doors. Press the "Tare" or "0/T" button to **zero the balance**. Zeroing the balance with the beaker or boat in place will make the balance **subtract out the flask or boat weight for you** when you add the liquid.
3. Open one of the doors and remove the flask or boat. Pour the liquid into the flask or boat, then carefully place it onto the center of the balance pan.
4. Close the glass door, wait for the balance reading to stabilize, then **write down all digits** displayed by the balance. Also **write down the units**, typically grams.
5. Remove your paper or boat with your liquid and shut the glass doors before leaving. If you have spilled any liquid on the balance, dry it off. Liquids can corrode the balance pan and destroy the balance if not removed!

If your liquid is volatile (evaporates easily), you may need to stopper the flask before weighing. Just make sure to zero the balance with the empty, stoppered flask.

A note on how precisely to measure chemicals: If you are told to measure 5.0 grams of powder, don't try to measure exactly 5.0000 grams. Just worry about matching the digits you are actually asked to measure. For example, a mass of 5.0493 g or a mass of 4.9983 grams (rounds to 5.0) would satisfy the directions. Write down **all the digits** from the balance when you make the measurement, since you may need all the digits for a calculation.

Measuring volume with a beaker

A beaker is essentially a cup, and is primarily designed for holding liquids rather than measuring them. Many beakers have printed scales on the side that can be used for "quick and dirty" measurements, and we will occasionally use them for this purpose in our lab.

When you pour a liquid into a container, you might observe that the liquid takes the shape of the container except on the top surface. The top surface **appears** flat. If you look closely, the liquid surface is actually slightly curved. You can most easily see the curvature at the the edge of the liquid surface. Liquids that are similar to water (almost everything you will work with in this course) curve up at the edges, forming a **meniscus**. When reading a volume measuring device, always read the scale at the **bottom of the meniscus**. By convention, we design scales on volume measuring devices to be read this way.




Illustration 2 - Cartoon diagram illustrating how to read the scale of a volume-measuring device

When reading a scale, always try to estimate a decimal place beyond the divisions on the scale. In the example in Illustration 2, we read the scale as 4.0 units, since the bottom of the meniscus appears to be right at the "4" mark. If it appeared to be between 4 and 5, we would try to guess how far it appeared to be and estimate the last digit. Numbers read from a scale should always include one estimated digit.

Measuring volume with a graduated cylinder

One of the most common volume measurement devices in the chemistry laboratory is the graduated cylinder - simply a glass cylinder printed with a scale.



Illustration 3 - Graduated cylinders in several colors and sizes

Cylinders have a scale that is more detailed than one on a beaker. Otherwise, they are read the same way.




Illustration 4 - Cartoon diagram of part of a graduated cylinder scale

In Illustration 4, we can see that the bottom of the meniscus lies between 1.7 and 1.8 units. It appears to be approximately halfway between these marks, so we record a volume of 1.75 units.

Equipment for the experiment

You will need a **50 mL graduated cylinder**, a **50 mL beaker**, a **50 mL Erlenmeyer flask**, a **#1 rubber stopper**, a **400 mL beaker**, a **Pasteur pipet with bulb**, and a **digital thermometer**.

Procedure

Important: Use the same analytical balance for the entire experiment!

Setting up

Fill the 400 mL beaker with deionized water. Allow the beaker to stand on your benchtop while you perform the next part of the experiment. Letting the water sit will allow it to come to a constant temperature.

The analytical balance

Weigh the 50 mL Erlenmeyer flask with the rubber stopper on the analytical balance **four times**. Make sure to re-zero the balance after each mass reading. Record the masses on page 9.

Density of water

Using a digital thermometer, measure the temperature of the water in the 400 mL beaker you set aside earlier. Record the Celsius temperature on the data sheets. Use water from this 400 mL beaker for the rest of the experiment.

The 50 mL beaker

Put 20 mL of deionized water into the 50 mL beaker. Use a Pasteur pipet to adjust the water level so that the bottom of the meniscus appears to lie on the 20 mL line. Go over to the analytical balance, and put your 50 mL **flask** with stopper in the center of the pan. Zero the balance - **with the flask inside** - then pour the water from the beaker into the flask. **Do not re-zero the balance at this point**. With the stopper in place and the glass doors shut, measure the mass of the water and record on page 9. **Empty the flask** into the sink and make sure the **outside** of the flask and stopper is dry. Repeat this entire section three times - with newly measured portions of water - for a total of four measurements of mass.

The flask should be stoppered during the zeroing of the balance and the actual weighing. The stopper will minimize any water loss due to evaporation during the measurement.

The 50 mL graduated cylinder

Put 20.0 mL of deionized water into the 50 mL graduated cylinder. Use a Pasteur pipet to adjust the water level so that the bottom of the meniscus appears to lie on the 20.0 mL line. Go over to the analytical balance, and put your 50 mL **flask** with stopper in the center of the pan. Zero the balance - **with the flask inside** - then pour the water from the cylinder into the flask. **Do not re-zero the balance at this point**. With the stopper in place and the glass doors shut, measure the mass of the water and record on page 9. **Empty the flask** into the sink and make sure the **outside** of the flask and stopper is dry. Repeat this entire section three times - with newly measured portions of water - for a total of four measurements of mass.

The flask should be stoppered during the zeroing of the balance and the actual weighing. The stopper will minimize any water loss due to evaporation during the measurement.

Data Analysis

First, look at the measurements of the mass of the empty stoppered flask from the analytical balance. These should vary by only about plus or minus 0.0001 grams. The balance manufacturer claims a precision of +/- 0.0001 grams, and you should be able to achieve this precision with simple mass measurements. (The analytical balance produces the most precise, or reproducible, measurements of any instrument in our laboratory!)

Using the mass of each water sample, you will calculate the volume of water dispensed. If you know

the density of water (its mass per unit volume), then it is a simple matter to change the measured mass into a volume. This is how you will determine if the beaker and cylinder dispense the claimed amount of water.

Calculate the volume of each portion of water dispensed by the cylinder and beaker from the masses you recorded using this formula.

$$\text{Volume} = \frac{\text{Mass}}{\text{Density}}$$

Look up the density of water in Table 1 for the temperature of water you recorded at the beginning of the experiment and use that for all your calculations. Write the results for these calculations on page 10. Provide a sample calculation (showing the numbers plugged into the formula you used) in the space provided..

Now calculate the average volume dispensed by the beaker and the average volume dispensed by the cylinder. Compare these averages, remembering that we were supposed to dispense exactly 20 mL in each case. Which measuring device gives a more accurate (closer to 20 mL) result?

Now look at the individual measurements. Is the variation between measurements larger for the beaker measurements or the cylinder measurements?: Which measuring device provides you with a more precise (reproducible) volume?

Make sure the names of all your group members are filled in on page 9 and that you answer all questions on page 11, then turn pages 9, 10, and 11 in to your instructor.

Table 1 - Density of water at different temperatures*.

<i>Temperature (°C)</i>	<i>Density (g/mL)</i>
16	0.99895
17	0.99878
18	0.99860
19	0.99841
20	0.99821
21	0.99800
22	0.99777
23	0.99754
24	0.99730
25	0.99705
26	0.99679
27	0.99652
28	0.99624
29	0.99595

*Adapted from data at http://www.ncsu.edu/chemistry/resource/H2Odensity_vp.html

Group information*Lab team information*

- Member 1: _____
- Member 2: _____
- Member 3: _____
- Date of experiment: _____

Data*The analytical balance*

Trial	Mass of dry stoppered 50 mL flask (g)
1	
2	
3	
4	

Density of water

Measured water temperature (C)	Density of water from table on page 8 (g/ml)

The 50 mL beaker

Trial	Mass of 50 mL portion of water, measured by beaker scale (g)
1	
2	
3	
4	

The 50 mL graduated cylinder

Trial	Mass of 50 mL portion of water, measured by cylinder scale (g)
1	
2	
3	
4	

Calculations

The 50 mL beaker

Fill in the chart below. Copy the masses from the data sheets, and calculate the volumes.

<i>Portion</i>	<i>Mass (g)</i>	<i>Calculated Volume (mL)</i>
1		
2		
3		
4		

- Average volume of portions from the beaker: _____

Sample calculation: Show in the space below how you calculated the volume of the first portion of water from the beaker.

The 50 mL graduated cylinder

Fill in the chart below. Copy the masses from the data sheets, and calculate the volumes.

<i>Portion</i>	<i>Mass (g)</i>	<i>Calculated Volume (mL)</i>
1		
2		
3		
4		

- Average volume of portions from the cylinder: _____

Sample calculation: Show in the space below how you calculated the volume of the first portion of water from the graduated cylinder.

